



ALA BIOLAB S.R.L.  
 www.alabiolab.ro  
 Bucharest, Street FETEȘTI, No. 54-56, Bl. 2,  
 postal code 032556  
 CUI 48680912  
 J40/15789/2023  
 E-mail : [ala@alabiolab.ro](mailto:ala@alabiolab.ro)  
 IBAN : RO64INGB0000999914337814

## 3D Rotating Mixer



LED display, the instrument is compact and convenient, and can be used in test benches, cold storage or incubators

- Three working modes can be switched: rotation mixing, shaking mixing and angle mixing
- Shell and panels are specially designed for easy cleaning and decontamination
- With multiple types of centrifuge tube clamps, which can be combined arbitrarily and easy to disassemble
- Through the adjustment of the rotating panel, the product can be rotated horizontally or vertically or at any angle
- With memory function, when turning on the instrument, it will display the last speed/angle setting value
- Widely used in biomolecular science, histochemistry, biochemistry, clinical medicine, etc.



<b>Model</b>	<b>TM40</b>
Speed range	10~40rpm
Movement mode	Rotation/Shake/Angle
Working mode	Continuous/Timed
Display	LED
Power	4W
Application tubes	0.5/1.5/2/5/7/10/15/50ml
Input voltage	100~240V 50/60Hz
Dimension	L240×W140×H210mm
Weight	1.2kg