



ALA BIOLAB S.R.L.
 www.alabiolab.ro
 Bucharest, Street FETEȘTI, No. 54-56, Bl 2,
 postal code 032556
 CUI 48680912
 J40/15789/2023
 E-mail : ala@alabiolab.ro
 IBAN : RO64INGB0000999914337814

Microplate Incubator Heating Shaker

| | | |
|---|--|---|
|  | <p>Low noise, elegant appearance, easy operation Continuous or timed operation with audible alarm and automatic cut-off Integrated bi-directional heating and orbital shaking provide maximum efficiency for incubation Provides three combination functions, including: Up to 4 microplates can be loaded at one time (96-well/384-well plate) Cell culture plate (24-well plate plate 48-well plate 96-well plate, etc.) Microplate Shaker/Integrated Bench top Incubator without Shaker Function Microplate Shaker Operating at Room Temperature Heated lid design with full top cover and flat heating plate provide bi-directional heating, ensures superior temperature uniformity and stability while also eliminating the effects of condensation</p> |  |
|---|--|---|

| | |
|------------------------------|-------------------------------------|
| Capacity | 4 pcs microplates or culture plates |
| Heating method | Heating film |
| Speed | 200~1500rpm |
| Speed accuracy | ±5% |
| Rotation amplitude | 2mm horizontal rotation |
| Time range | 1~99h59min |
| Temperature control range | 1~99h59min |
| Temperature control range | Room temperature +5C~70°C |
| Temperature uniformity | ±0.5% |
| Temperature dispaly accuracy | 0.1°C |
| Heating time | ≤10min |
| Power | 250W |
| Weight | 8kg |
| Input voltage | 220V 50Hz |
| Dimension | L340×W210×H130mm |